

Let Clarity bring your world into focus.



Compatible with desk and laptop computers.  
Includes remote control.

**CLARITY**

1-800-575-1456  
www.clarityusa.com

## ClarityDeskmate

For the person on the go, the Clarity *Deskmate* system offers a high-quality mobile video magnification solution. Deskmate lets you travel easily to work or school with minimal equipment to carry.

2-55X<sup>MAGNIFICATION</sup>

*Deskmate* and *Deskmate Plus* combines Clarity's Flex video magnifier with a 17" high-resolution LCD monitor for excellent visibility. *Deskmate* is ready to go wherever you go. Setup is a snap — just plug in and switch on!

*Deskmate* features Clarity's unique Vector Imaging System™ for smooth, hands-free automatic focusing. You'll view three-dimensional items with amazing sharpness, and see your favorite photos in full color!

Our versatile rotating camera design provides crystal-clear, full color viewing at any range: close-up, intermediate or far away. Adjusting from near to distant viewing is quick and simple with *Deskmate's* tilt feature. With its convenient remote control, using *Deskmate* couldn't be easier.

*Deskmate Plus* integrates a rechargeable battery in the base of the unit which gives the user additional independence.

Clarity's **LIFETIME WARRANTY** gives you peace of mind while you're on the move. Like all Clarity products, *Deskmate* is backed by our superior customer service, and is individually handcrafted in the USA.

Let **Clarity** bring **the world into focus.**

Please note: The actual product may differ slightly in appearance from that shown.

**The most important thing to Clarity are our customers. Let's hear what they have to say ...**

"I am writing to express my enthusiasm for my new Clarity system. There are a hundred things I can do with my video magnifier. I can read the newspaper, books, and articles. I can magnify items with print I need to read, such as a vitamin bottle. I am able to do my bookkeeping, balance my check book, figure my taxes, write checks and look at pictures. ...Also, I have found that as a company, Clarity will do everything they say they will do and more."

Glenn Christy, low vision customer

## Lynx



For the ultimate in portability go to [www.clarityusa.com](http://www.clarityusa.com) or call 800-575-1456 to see the Lynx.

## Specifications

Magnification:

- 2x-55x

Features:

- Color, Black and White and Inverse
- Foldable and portable
- Mirror Mode
- Zoom, contrast, and mode controls
- Low battery indicator (Plus Model)

Unit Dimensions:

- Height = 14.75"
- Width = 17"
- Depth = 6.75" @ Base

Weight:

- 11lbs. (Deskmate)
- 14lbs. (Deskmate Plus)

LCD:

- 17" High Contrast
- 19" LCD (Deskmate Option Only)

Power Requirements:

- AC Adapter

Accessories

- AC adapter
- Carrying Case
- Remote control

**CLARITY**

**800-575-1456**

for more details.