

Let Clarity bring your world into focus.

NEW!



CLARITY

1-800-575-1456
www.clarityusa.com

ClarityLynx

The Lynx portable magnifier offers a compact design that will fit easily into any work space from a desktop to a countertop to fit any lifestyle. The self contained system offers 5 hours of battery life with over 10 inches of viewing space and weighs only 8 pounds.

3-36X MAGNIFICATION

Lynx is the latest in portable technology. It is a portable system with battery power for the ultimate in traveling flexibility, with distance viewing, and two camera options. Now you can go to the theatre, restaurants, business meetings, or classroom and have everything you need at your fingertips.

Lightweight, large screen, and long battery life makes this product a truly ideal solution for those of you who need more magnification than handheld magnifiers can offer. Choose color, positive or negative text, or yellow on blue viewing mode and adjust brightness to best suit your needs for optimal reading and writing results.

The handy mouse camera gives you total control for reading at a pace comfortable for you. The main camera for the unit gives you the ability to view pictures, objects, read and write, and even distance view with ease. The adjustable arm is made for you to set the camera to fit your personal viewing needs perfectly. A quick turn of the camera and you are distance viewing instantly. It even offers our Quick release option allowing you the freedom to create additional magnifiers to fit your lifestyle.

Please note: The actual product may differ slightly in appearance from that shown.

The most important thing to Clarity are our customers. Let's hear what they have to say ...

"I'm a teacher of the Visually Impaired Program who recently helped a junior high student obtain a Clarity. Without having had a demonstration or any prior experience with the Clarity, I was impressed with the ease in which I quickly hooked up the device and had it operational.

Likewise, when instructing the student on the device and process for using the device and moving it from setting to setting - the student was able to efficiently perform the process within one session."

Kathy Oliver

Financing is available.
Call 800-575-1456 to find out more details.

Specifications

Magnification:

- 3x-36x

Features:

- Color, Black and white, Inverse, Yellow on Blue, and black on yellow
- Foldable and portable
- Self-viewing mode
- Power save
- Read, write, distance view
- Quickrelease camera with autofocus
- USB output
- Mouse camera
- Zoom, contrast, and mode controls
- Low battery indicator

Unit Dimensions:

- 14" x 9"

Weight:

- 8 lbs.

LCD:

- 10.4" High Contrast

Power Requirements:

- Five hour battery life

Accessories

- AC adapter
- Carrying Case
- Mouse camera
- USB output cable
- Remote control