

Flexmate

The Flexmate video magnifier is expertly designed for reading, writing, and distance viewing. This makes the Flexmate a perfect solution for students, working professionals, hobbyists or homemakers. Clarity makes it possible for you to add options to fit your lifestyle!

2-60X MAGNIFICATION



Shown with optional XY Table

The Flexmate offers:

- 14" adjustable flex arm with clamp mount*
- Quick Release Camera with automatic focus
- Read, write and distance viewing
- Tilt feature for self-viewing
- Tactile one touch zoom and contrast
- Color, black and white, or inverse
- Focus lock
- Clamp and Carry bag for portability
- Remote control

Advanced Features:

- Masking • Manual Focus • Indoor and Outdoor Lighting Options • Unique Media Viewing Ability



"I have Macular Degeneration. I'm a portrait painter, but I had to give up painting. Now with my Clarity Flex I'm able to enlarge a photograph or a live model and paint away!"

* Optional Long Arm Flex with adjustable 30" arm.

Carrymate

The Carrymate is the ultimate portable video magnifier and is ideally designed to fit easily into any work space. It's fully self-contained with a rechargeable battery. The Carrymate has no additional cables so it's ready to go with you anywhere.

4-36X MAGNIFICATION



The Carrymate offers:

- 10.4" LCD screen
- Base footprint 10"x12"
- 9 hour battery with charge indicator
- Expandable/collapsible arm
- Quick Release Camera with automatic focus
- Read, write or distance viewing
- Tilt feature for self-viewing
- Tactile one touch zoom and contrast
- Color, black and white, or inverse
- Focus Lock
- Hard carry case for portability
- Remote control

Advanced Features:

- Masking • Manual Focus • Indoor and Outdoor Lighting Options • Unique Media Viewing Ability



"I use the Carrymate in a classroom and library. It's a great product. It allows me to see the black board and my reading speed is better now. I definitely recommend the Carrymate!"

PCMate

If you're a Laptop user the PCMate is for you! USB connection to your PC means no batteries, cords, or power outlet needed. Included Capture software offers split screen, full screen, and capture ability to view images later.

4-60X MAGNIFICATION




Product does not include laptop.

The PCMate offers:

- 9" adjustable flex arm
- USB interface for PC compatibility
- Capture images, split screen, window sizing
- Camera with automatic focus
- Read, write and distance viewing
- Tilt feature for self-viewing
- Tactile one touch zoom and contrast
- Color, multiple color options, black and white, or inverse
- Focus lock
- Carry bag for portability
- Remote control

Advanced Features:

- Masking • Manual Focus • Indoor and Outdoor Lighting Options • Unique Media Viewing Ability



"Before the PCMate, I strained hard to view the speaker, marker-board, and power point slides. I feel so much more productive and involved in my corporate meetings now. Best of all, the PCMate can be easily adjusted to read the hand outs given at meetings!"

Junior

This state of the art hand-held video magnifier can be taken to your favorite restaurant, the grocery store or anywhere else your busy life takes you. You will have instant video magnification that is extremely portable and can fit easily into a purse or backpack.

3-9X MAGNIFICATION

The Junior offers:

- 4" LCD screen
- 11 ounces, 5 x 5 x 1.5 inches
- 3 hour battery life
- Reading and writing ability
- Standard wall charger*
- Video output to use with larger screen
- Contrast and Mode controls
- Color, Black and White, Inverse, and Color Select mode
- Carry case for portability

* Optional car charger



"The Junior makes it so easy to read menus at restaurants and price tags when I'm shopping!"



"There are a hundred things I do with my Clarity Discovery every day! I can read the newspaper, books, trade journals, prescription pill containers, balance my checkbook and look at photographs. I heartily recommend Clarity!"

Discovery

If your primary needs for magnification include reading, writing, and other close-up tasks, the Clarity Discovery is the video magnification solution for you.

4-60X MAGNIFICATION

The Discovery offers:

- Easy to use controls
- Adjustable magnification
- One touch power switch
- Full color, black and white, or inverse options
- Contrast control
- Autofocus
- High contrast 17" LCD
- Full size XY reading table
- Monitor tilt feature
- Read, write or view pictures
- Black or white base




Lifetime Warranty

Lifetime warranty non-transferable. One year warranty for batteries.



2222 Park Place, Suite 1C
Minden, NV 89423
Web: www.clarityusa.com

Tel: 800-575-1456
Fax: 775-783-0966
Email: info@clarityusa.com



"The fifth grade student was ecstatic! She could see her teacher's face, the board, the overhead and her classmates for the first time! Thank you for creating Deskmate."

Deskmate

Deskmate combines Clarity's Flex video magnification with a 17" high-resolution LCD monitor. Lightweight and portable, the Deskmate is ready to go when you are. Quick and easy setup, just plug it in and switch on.

2-55X MAGNIFICATION

The Deskmate offers:

- 17" LCD screen
- Collapsible for easy transport
- Expandable/collapsible arm
- Quick Release Camera with automatic focus
- Read, write and distance viewing
- Tilt feature for self viewing
- Tactile one touch zoom and contrast
- Color, black and white, or inverse
- Focus lock
- Carry bag for portability
- Remote control

Advanced Features:

- Masking • Manual Focus • Indoor and Outdoor Lighting Options • Unique Media Viewing Ability

WWW.CLARITYUSA.COM

800-575-1456

WWW.CLARITYUSA.COM