

# *Press Information*

FOR IMMEDIATE RELEASE

## Industry's Most Portable and Complete CCTV for Classroom and Executive Use

### Note to Editor:

- ♦ Lynx is a highly portable digital magnifier designed for classroom and business use.
- ♦ Lynx allows the user flexibility between desktop and distance viewing.
- ♦ Lynx provides magnification between 3x - 36x.
- ♦ Lynx folds up like a laptop, weighs only 8 pounds and can be taken between classes or meetings.
- ♦ Lynx includes the unique MouseCam for quick desktop reference while in distance mode.



LIVERMORE, CA - February 3, 2008 - Clarity introduces Lynx, the industry's most portable and complete CCTV for classroom and executive use. The Lynx offers a compact design that will fit easily into any work space from a desktop to a countertop to fit any lifestyle. The self contained system offers 5 hours of battery life, over 10 inches in viewing space and weighs only 8 pounds. According to Dave Brody, CEO of Clarity, "The Lynx gives the student or business executive the tools they need to be successful in a highly portable, compact design. The fact that it is the size and weight of a small laptop makes it ideal."

### Like a Laptop

Traditionally video magnifiers have been larger and cumbersome to use in school or business. When students need to move from classroom to classroom the need for higher portability becomes essential. These same requirements apply to the business executive. The Lynx provides both the opportunity to be more effective in their day to day activities without the added burden of carrying a bulky video magnifier.

### Summary of Key Features

- ♦ Weighs only 8 lbs. and folds up into laptop type carry case for easy portability.
- ♦ Provides magnification between 3x-36x.
- ♦ Quick release camera with autofocus.
- ♦ Ability to toggle between main camera and MouseCam.
- ♦ USB connection for image capture via PC.
- ♦ Color, B&W, Inverse, Yellow on Blue and Black on Yellow modes allow for easier viewing.
- ♦ Power save/Auto shut off maximize the 5 hours of battery life.
- ♦ Low battery indicator lets you know when it is time to recharge.



Editors Contact: [info@clarityusa.com](mailto:info@clarityusa.com)

Customer Service: 1-800-575-1456