

# Comoframe



## The Trial Frame

Refraction and the following control of visual acuity can be simple when the person who undergoes the test is at ease and collaborates. On the contrary it becomes difficult if the environment, the lights and the situation cause apprehension and intolerance.

The problem increases with babies. Dynamic and unpredictable they become suspicious if any device prevents the peripheral vision.

They get nervous any try look sideways to check what is being done.

Hence the need to use simple devices for the various visual tests.

The trial frame COMOFRAME is the result of experiences carried out by ophthalmologists, strabologists, optometrists and orthoptists.

Light and comfortable they are accepted by children as a normal pair of glasses. It is always well balanced on the face because it is held at the back by elastic and allows the tester to have the patient's full concentration.

The lens holder rings use standard 38 mm trial lens stes and can be supplied, on request, with Tabo or international indications.

The trial frame COMOFRAME is available in a small handy hard case.



### Available colors:

<input type="checkbox"/>	01 - Bianco Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	02 - Rosa Pallido Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	03 - Azzurro Pallido Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	04 - Giallo Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	05 - Rosa Trasparente	<input type="checkbox"/>
<input checked="" type="checkbox"/>	06 - Azzurro Trasparente	<input checked="" type="checkbox"/>
<input type="checkbox"/>	07 - Verde Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	08 - Blu Trasparente	<input type="checkbox"/>
<input type="checkbox"/>	14 - Rosso	<input type="checkbox"/>
<input type="checkbox"/>	16 - Verde Scuro	<input type="checkbox"/>
<input type="checkbox"/>	17 - Grafite	<input type="checkbox"/>
<input type="checkbox"/>	18 - Grigio Perla	<input type="checkbox"/>
<input type="checkbox"/>	19 - Amaranto	<input type="checkbox"/>
<input type="checkbox"/>	20 - Marrone	<input type="checkbox"/>

**Patent:** N. 21555 B/85

### Accessories:

- special elastic
- handy hard case

