

The 23g Transconjunctival Vitrectomy System from Eye Technology

- 3 x Universal Cannulas
- 2 x Closure Valves
- 3 x Blunt Insertors
- 1 x Universal Infusion Line
- Eye Tech Cannula Forcep
- Pressure Plate
- Sterilisation Box

The new Transconjunctival 23g Vitrectomy System from Eye Technology marks another milestone in product design. The system incorporates a unique self closing valve system making it possible to perform transconjunctival sutureless vitreous surgery easier and more efficiently than ever before. Good surgical indications such as epiretinal membrane proliferation, macular oedema, macular hole surgery, and simple vitreous haemorrhage, are ideal for the 23-gauge vitrectomy system. Together with the new range of 23g instruments available from Eye Technology, more complex cases, such as rhegmatogenous retinal detachment, proliferative vitreoretinopathy, and diabetic retinopathy can be treated more easily.



Eye Technology
www.eye-tech.co.uk

19 Totman Crescent, Rayleigh SS6 7UY, UK.
Tel: +44 (0)1268 771949 Fax: +44 (0)1268 771349
Email: sales@eye-tech.co.uk



Part No: 50-0101
23g Conjunctival Trocar Set


23 GAUGE ACCESSORIES




Part No: 50-0102
23g Forceps




Part No: 50-0103
Blunt Insertor



Part No: 50-0106
Closure Valve



Part No: 50-0104
PVC Infusion Line



Part No: 50-0107
Silicon Infusion Line



Part No: 50-0105
Pressure Plate



Part No: 50-0108
23g MVR Knife