


MEDICOSPOT

MEDICOSPOT MS

- Produces brilliant and clear spot light of .3cm diameter at the distance of 15cm.
- Size of spot light variable by rotating ring.
- Used with 2pcs. of batteries.
- Contents: MS, Spare bulb.
- Bulb: MS bulb (L-12)
- Instrument Weight: 90grm
- Total Weight: 140grm
- Case Size: 180mm(D) × 90mm(W) × 50mm(H)



MS

TRANSILLUMINATOR MS II

- Pocket type transilluminator.
- Used for transillumination in such fields as Ophthalmology, obstetrics, gynecology and surgery.
- Used with 2ps. of batteries
- Contents: MS II, Spare bulb.
- Bulb: Laryngoscope bulb (L-06)
- Instrument Weight: 90grm
- Total Weight: 140grm
- Case Size: 180mm(D) × 90mm(W) × 50mm(H)



DOCTOR LIGHT DL

- Slim and feather-weight pocket light.
- Sharp and brilliant spot light produced by a two-element condensing lens system.
- Contents: DL, spare bulb.
- Bulb: DL bulb (L-37)
- Instrument Weight: 50grm
- Length: 130mm
- Diameter: 14mm
- Case Size: 170mm(D) × 60mm(W) × 20mm(H)



DL